# Duroval® C

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
## Reagent to determine the Carbonate Hardness in water

The Duroval-C offers a simple and rapid method to determine Carbonate Hardness, based on the principle of colour change.

Each kit contains:

- 1 calibrated test tube
- 1 stopper to close the test tube
- 1 dropper bottle C-indicator
- 1 bottle containing titration solution
- 1 test pipette, graduated in mmol/l and °dH.

The C-indicator disperses in hard water to blue-green. By adding the titration solution, the carbonate hardness is neutralized, and the colour turns orange. The amount of titration solution used measures the carbonate hardness contents of the water sample. The carbonate hardness value can be read off on the graduated test pipette.

### Directions:

1. Rinse out the test tube with the water sample and fill up to the 5 ml mark (fig. 1).
2. Add 1 drop of C-indicator and disperse by swirling (fig. 2). If the sample turns to orange, no carbonate hardness is present.
3. If the sample turns to blue-green, close the test tube with the stopper (fig. 3).
4. Aspirate titration solution from the storage bottle with the graduated pipette (fig. 4 a) until the lower rim of the black piston ring is at the same level as the 0-mark (uppermost mark on the graduated scale) (fig. 4 b). For an accurate test result, eliminate any air bubbles by depressing the plunger fully with its contents into the storage bottle, keeping the tip of the pipette immersed in the solution. Repeat the process if necessary to evacuate all air bubbles. Finally, refill the pipette to the graduated 0-mark with the titration solution.
5. With a slight twist firmly insert the pipette tip into the orifice of the test tube stopper (fig. 5 a). Whilst shaking (fig. 5 b), depress the pipette plunger very slowly (fig. 5 c) until the colour of the water sample turns from blue-green to orange (fig. 6 a).
6. Read the degree of carbonate hardness taken at the lower rim of the black piston ring on the graduated pipette (fig. 6 b) and take note. If the contents of one pipette are insufficient to effect the colour change, refill the pipette according to 4 with fresh titration solution up to the 0-mark as described and continue the titration in the same water sample according to 5 and 6.  
Note! When determining the carbonate hardness, allow for the quantity of titration solution previously used.

$$1 \text{ } ^\circ\text{dH} = 1,79 \text{ } ^\circ\text{f} = 17,9 \text{ mg/l CaCO}_3 = 0,36 \text{ mmol/l}$$

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