testoval[®]

KUPFER-TEST

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Directions for Use

Remove the test cell from the comparator and fill it up to the mark with water sample. Add a dosing spoon of copper - reagent and dissolve by stirring with the red spoon. 2 minutes after the reagent is completely dissolved, reinsert the test cell into the comparator, hold the comparator against a bright background or a light source, and match the colour intensity of the central field (test solution) with one of the colours in the peripheral colour fields. Read the value in ppm Cu, which is printed under the matching colour field.

Should the colour in the central field be more intensive than any of the colours in the peripheral fields, repeat the test with a water sample diluted 1 : 10. For this purpose fill the graduated plastic tube up to the lower mark with sample water and dilute with copper-free water (e. g. distilled water) to the upper mark. Transfer the diluted sample into the test cell, which has been rinsed with copper-free water and proceed as described above. The value obtained must in this case be multiplied by 10.

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